



## Total Sulfite Assay Kit (Colorimetric)

### Product Information

---

**Cat**

Kit-1025

---

**Common Name**

Sulfite

---

**Cat.No.**

Kit-1025

---

**Description**

Sulfites/sulphites are substances that occur naturally in the human body due to the metabolism of amino acids containing sulfur in their side chains. They are readily oxidized to sulfates via enzymatic reactions and excreted in urine at the rate of 1000 mg/day. Sulfites are often considered as allergens due to adverse symptoms observed in asthmatic patients. Exogenous sources of sulphites include polluted air, food and beverages containing sulphur dioxide (SO<sub>2</sub>). SO<sub>2</sub> is a molecule that easily reacts with several small compounds including aldehydes, ketones, anthocyanins, cobalamine, thiamine, NAD, flavins, among others. SO<sub>2</sub> also interacts with cysteine residues in proteins and promotes cross-linking. Sulfites are also used as regulated food preservatives/additives/enhancers in dried fruits, wine, beer, etc. They are considered as GRAS (Generally Recognized as Safe) in food and beverages. In wine, SO<sub>2</sub> reacts with sugars, aldehydes and anthocyanins. The term sulfite is used for all sulfite-derived molecules: bisulfites, metasulfites, and SO<sub>2</sub>. Total sulfite is defined as the sum of bound and free sulfite. Total Sulfite Assay kit is a simple and sensitive assay to detect small concentrations of sulfites in a variety of food and beverage samples. This assay is based on the oxidation of sulfite to sulfate producing a stable signal at 570 nm, which is directly proportional to the amount of sulfite in the sample. This assay is very sensitive and can detect as low as 20 μM of sulfite in a variety of samples.

---

**Applications**

Estimation of sulfite in various food and beverage samples.

---

**Storage**

---



## Total Sulfite Assay Kit (Colorimetric)

-20°C

---

### Shipping

Gel Pack

---

### Size

100 assays

---

### Kit Components

Sulphite Assay Buffer; Sulphite Stabilizer; Sulphite Probe; Sulphite Oxidizing Mix; Sulphite Enzyme Mix; Sulphite Standard

---

**Detection method** Absorbance (570 nm)

---

### Features & Benefits

Highly sensitive, quick, colorimetric assay.

Can detect as low as 20 µM of sulfite in a variety of samples.

---