

CECR2 Inhibitor Screening Assay Kit

Product Information

Cat.No.

Kit-1707

Size

384 reactions

Description

The CECR2 Inhibitor Screening Assay Kit is designed to measure the inhibition of CECR2 binding to its substrate. The key to the CECR2 Inhibitor Screening Assay Kit is the highly specific binding of the CECR2 bromodomain to the acetylated histone substrate. With this kit, only three simple steps on a microtiter plate are required. First, a sample containing CECR2 and an inhibitor of choice is incubated with the biotinylated substrate for thirty minutes. Next, acceptor beads are added, then donor beads, followed by reading the Alpha-counts.

Applications

Useful for the study of bromodomain binding assays, screening bromodomain inhibitors and selectivity profiling.

Storage

Stable at least 12 months from date of receipt, when stored as directed. Kit components require different storage conditions. Be sure to store each component at the proper temperature upon arrival.

Kit Components

His-CECR2: 100 µg; -80°C BET Bromodomain Ligand: 400 µl; -80°C Non-acetylated Ligand 1: 200 µl; -80°C 3x BRD Homogeneous Assay Buffer 2: 4 ml; -20°C 3x BRD Homogeneous Detection Buffer 2: 3 ml; -20°C