



Live-Dead Cell Staining Kit

Product Information

Cat

Kit-2360

Cat.No.

Kit-2360

Product Overview

The assay range is 0.156–10 ng/ml ANGPTL3/ml. The lowest level of ANGPTL3 that can be detected by this assay is 75 pg/ml.

Applications

Stained live and dead cells can easily be visualized by fluorescence microscopy.

Storage

-20°C

Shipping

Gel Pack

Size

100 stainings

Kit Components

Solution A (1 mM Live-Dye); Solution B (1 mg/ml PI); Staining Buffer

Target Species

Mammalian

Detection method Fluorescent microscopy (Ex/Em488/518 nm) to detect stained live cells; (Ex/Em 488/615) to detect stained dead cells

Features & Benefits

Simple procedure; takes less than 20 minutes; Fast and convenient