



Fluorescein Active Caspase-9 Staining Kit

Product Information

Cat

Kit-2385

Cat.No.

Kit-2385

Product Overview

Activation of caspases plays a central role in apoptosis. The CaspGLOW;Fluorescein Active Caspase-9 Staining Kit provides a convenient means for sensitive detection of activated caspase-9 in living cells. The assay utilizes the caspase-9 inhibitor LEHD-FMK conjugated to FITC (FITC-LEHD-FMK) as a fluorescent marker. FITC-LEHD-FMK is cell permeable, nontoxic, and irreversibly binds to activated caspase-9 in apoptotic cells.

Applications

Sensitive detection of activated caspase-9 in living cells.

Storage

-20°C

Shipping

Gel Pack

Kit Components

FITC-ATAD-FMK; Wash Buffer; Z-VAD-FMK

Target Species

Mammalian

Detection method Fluorescence microscopy, Flow cytometry and fluorescence plate reader (Ex. = 485 nm and Em. = 535 nm)

Features & Benefits

Simple one-step procedure; takes only 1-2 hours; Fast and convenient; The FITC label allows for



CREATIVE **BIOMART**[®]
Assay Kit

Fluorescein Active Caspase-9 Staining Kit

direct detection of the activated caspases in apoptotic cells

Tel: 1-631-559-9269 1-516-512-3133

Fax: 1-631-938-8127

Email: info@creative-biomart.org

45-1 Ramsey Road, Shirley, NY 11967, USA