



## Citric Acid (OA) Content Assay Kit

### Product Information

---

#### Product Overview

---

Oxalic acid is a dicarboxylic acid commonly found in plant tissues. As an antagonist of mineral elements, consuming vegetables with high oxalic acid content over a long period can easily lead to diseases such as arthritis, hypocalcemia, and kidney stones. Oxalic acid can serve as a complexing agent, bleaching agent, precipitating agent, and reducing agent, and it is widely used in industries such as pharmaceuticals, dyeing, and plastics. Fe<sup>3+</sup> can react with sulfosalicylic acid under acidic conditions to form a purple complex, which has a characteristic absorption peak at 510 nm. Oxalic acid and oxalate can lighten the color of this purple complex, and the absorbance value of the Fe<sup>3+</sup> and sulfosalicylic acid complex decreases with the increase of oxalic acid. By detecting changes in absorbance, the oxalic acid content can be determined.

---

#### Size

---

48 Samples

---

#### Storage

---

4□

---

#### Shipping

---

Ice pack

---

#### Kit Components

---

Extraction Liquid: Liquid, 60mL × 1 bottle, store at 4°C.

Reagent One: Powder, mg × 1 vial, store at 4°C. Before use, shake to bring the powder to the bottom, then add 1.5mL of distilled water to dissolve and reserve.

Reagent Two: Liquid, 20mL × 1 bottle, store at 4°C.

Reagent Three: Liquid, 1.5mL × 1 vial, store at 4°C.



CREATIVE **BIOMART**<sup>®</sup>  
Assay Kit

## Citric Acid (OA) Content Assay Kit

Standard: Powder, mg × 1 vial, store at 4°C. This reagent is used if recalibration of the standard curve is necessary.

---

### **Materials Required but Not Supplied**

---

Microplate reader, 96-well plate, benchtop centrifuge, adjustable pipette, water bath/metal bath, mortar, ice, distilled water.

---

Tel: 1-631-559-9269 1-516-512-3133

Email: [info@creative-biomart.org](mailto:info@creative-biomart.org)

Fax: 1-631-938-8127

45-1 Ramsey Road, Shirley, NY 11967, USA