



Furantoin Metabolite Detection Kit

Product Information

Product Overview

Drug / Contaminant Residue Detection Kits

Description

Detects AMOZ, a metabolite of Furantoin, in food samples such as meat, milk, and aquatic products.

Size

48t/96t

Applications

- Monitoring of AMOZ residues in animal-derived foods
 - Regulatory enforcement testing
 - Ensuring food export safety compliance
-

Features & Benefits

- Highly specific to Furantoin metabolites (AMOZ)
 - Rapid and easy-to-perform ELISA format
 - Suitable for meat, fish, and dairy products
 - Sensitive to low-level contamination
-