



Furazolidone Metabolite Detection Kit

Product Information

Product Overview

Drug / Contaminant Residue Detection Kits

Description

Detects AOZ, the key metabolite of Furazolidone, in meat, milk, and seafood products.

Size

48t/96t

Applications

- Detection of Furazolidone residues in food
 - Compliance testing for exports
 - Veterinary drug surveillance
-

Features & Benefits

- Fast and efficient AOZ detection
 - Suitable for various matrices including muscle and liver
 - ELISA format enables high-throughput testing
 - Supports residue monitoring programs
-