



## Aspartate Assay Kit

### Product Information

#### Cat.No.

Kit-0109

#### Product Overview

Aspartate can be quantified in the range between 0.1-10 nmoles/well (2-200  $\mu$ M). (100 assays)  
Simple procedure; takes less than ~40 minutes.Fast and convenient

#### Size

100 assays

#### Applications

Aspartate can be quantified in the range between 0.110 nmoles/well (2-200  $\mu$ M).

#### Target Species

Mammalian

#### Usage

This product is furnished for LABORATORY RESEARCH USE ONLY. Not for diagnostic or therapeutic use.

#### Kit Components

Aspartate Assay BufferProbe (DMSO solution)Serum CleanUp MixAspartate Enzyme MixConversion  
MixAspartate Standard (100 mM)

**Detection method** ColorimetricFluorescence Plate Reader

#### Compatible Sample Types

Cell Culture SupernatantsPlasmaSerumUrine