



## Fructose-6-Phosphate Assay Kit

### Product Information

#### Cat.No.

Kit-0331

#### Product Overview

The PicoProbe F6P Assay Kit can detect F6P in the range of 0.01 to 0.5 nmoles with detection sensitivity ~1  $\mu$ M of F6P. (100 assays) Simple procedure; takes ~ 35-40 minutes to complete the assayFast and convenientSensitive, efficient and accurate assay

#### Size

100 assays

#### Applications

The PicoProbe F6P Assay Kit can detect F6P in the range of 0.01 to 0.5 nmoles with detection sensitivity ~1  $\mu$ M of F6P.

#### Target Species

Mammalian

#### Usage

This product is furnished for LABORATORY RESEARCH USE ONLY. Not for diagnostic or therapeutic use.

#### Kit Components

F6P Assay BufferF6P ProbeF6P Enzyme MixF6P ConverterF6P DeveloperF6PStandard (10  $\mu$ mol)

**Detection method** Fluorescence Plate ReaderFluorometer

#### Compatible Sample Types

Cell Culture SupernatantsFood ProductsGrowth MediaPlasmaSerumUrine