



# Phosphatidylserine Apoptosis Kit (Orange Fluorescence)

## Product Information

### Cat.No.

Kit-0747

### Product Overview

Phosphatidylserine Apoptosis Kit (Orange Fluorescence) is a fluorescent assay used to monitor cell apoptosis.

### Description

There are a variety of parameters that can be used for monitoring cell viability. This particular kit is designed to monitor cell apoptosis by measuring the translocation of phosphatidylserine (PS). In apoptosis, PS is transferred to the outer leaflet of the plasma membrane. The appearance of phosphatidylserine on the cell surface is a universal indicator of the initial/intermediate stages of cell apoptosis and can be detected before morphological changes can be observed.

### Applications

This kit is a set of tool for monitoring cellular functions.

### Usage

For research use only (RUO)

### Storage

Keep at 4 °C and protect from light.

### Kit Components

Component A: Apopxin Orange (100X Stock Solution) 100  $\mu$ L  
Component B: Assay Buffer 10 mL

**Detection method** Fluorometric

### Compatible Sample Types

Adherent cells, Non-adherent cells

### Features & Benefits

Non-Radioactive: No special requirements for waste treatment. Continuous: Easily adapted to



CREATIVE **BIOMART**<sup>®</sup>  
Assay Kit

## Phosphatidylserine Apoptosis Kit (Orange Fluorescence)

automation with minimal hands on time. Convenient: All essential assay components are included. Optimized Performance: Optimized for detecting the translocation of phosphatidylserine (PS). Enhanced value: Less expensive than the sum of individual components.

---

Tel: 1-631-559-9269 1-516-512-3133

Email: [info@creative-biomart.org](mailto:info@creative-biomart.org)

Fax: 1-631-938-8127

45-1 Ramsey Road, Shirley, NY 11967, USA