



Annexin V-PE-Cy5 Apoptosis Detection Kit

Product Information

Cat

Kit-2389

Cat.No.

Kit-2389

Product Overview

The Annexin V-PE-Cy5 Apoptosis Detection Kit is based on the observation that soon after initiating apoptosis, most cell types translocate the membrane phospholipid phosphatidylserine (PS) from the inner face of the plasma membrane to the cell surface. Once on the cell surface, PS can be easily detected by staining with a fluorescent conjugate of Annexin V, a protein that has a strong natural affinity for PS. The one-step staining procedure takes only 10 minutes. In addition, the assay can be directly performed on; Live cells, without the need for fixation. Results can be analyzed by flow cytometry or fluorescence microscopy.

Applications

Detect early/middle stages of apoptosis; differentiate apoptosis from necrosis.

Storage

4°C

Shipping

Gel Pack

Kit Components

Annexin V-PE-Cy5; 1X Binding Buffer

Target Species

Mammalian

Detection method Flow cytometry (Ex = 488 nm; Em = 670 nm) and fluorescence microscopy

Features & Benefits



CREATIVE **BIOMART**[®]
Assay Kit

Annexin V-PE-Cy5 Apoptosis Detection Kit

Simple one-step procedure; takes only 10 minutes; Fast and convenient, no need of fixation of cells

Tel: 1-631-559-9269 1-516-512-3133

Fax: 1-631-938-8127

Email: info@creative-biomart.org

45-1 Ramsey Road, Shirley, NY 11967, USA